

The Necessity of Bleeding Control Kits in Schools

Include Bleeding Control Solutions in Your School Budget.

The right tools in place make a difference when minutes matter. Each kit is equipped with life-saving components.

100K 
Students attended schools where school violence incidents happened
**Source: SIEPR 2018-2019

140 
People affected by injuries or loss of life in 2022
**Source: Edweek 2022

5 Min. 
All the time it takes for a person to die from blood loss.**
**Source: facts.org, American College of Surgeons

Bleeding Control Kit Limb, Chest/Lung, & Torso Wounds



QTY.	DESCRIPTION
1	Bandage Shears
1	Bleeding Control Instructions
2	Vented Chest Seal
1	Emergency Blanket, 52" x 84"
1	Hemostatic Gauze Pad, 4" x 4"
1	Israeli Bandage
1	Marker
4	Nitrile Exam Gloves
1	Tourniquet
1	Trauma Pads, 5" x 9"
1	Windlass Tourniquet



Part No.: 91137

Bleeding Control Kit Limb & Torso Wounds



QTY.	DESCRIPTION
1	Windlass Tourniquet
1	Hemostatic Gauze Pad, 4" x 4"
1	Black Marker
2	Nitrile Exam Gloves
2	Trauma Pads, 5" x 9"
1	Bleeding Control Instructions



Part No.: 91483

Bleeding Control Kit Limb & Chest/Lung Wounds



QTY.	DESCRIPTION
1	Windlass Tourniquet
2	Vented Chest Seals
1	Black Marker
2	Nitrile Exam Gloves
2	Trauma Pads, 5" x 9"
1	Bleeding Control Instructions



Part No.: 91481

Discover Essential Life-Saving Tools



Windlass Tourniquet
Part No.: 91160
Quickly restrict blood flow on limbs and extremities
Treats limb wounds

A tourniquet is used to stop severe bleeding on an arm or leg when direct pressure alone does not. A tourniquet is wrapped around the limb above the wound and tightened to restrict blood flow.

Some tourniquets like the First Aid Only QuickTight Tourniquet have a built-in windlass wand that is twisted to cut off blood flow quickly and then clipped into a clasp.



Vented Chest Seal
Part No.: 91194
Preserves lung function for critical puncture wounds
Treats chest/lung wounds

A chest seal covers a lung puncture wound and allows the lung to continue to function. Vented chest seals have channels to allow excess air trapped in the chest cavity to vent out while also preventing air from the outside from being sucked into the chest cavity. The chest seal is coated with a thick gel to adhere to a victim's chest without interference from hair or moisture.



Hemostatic Gauze Pad
Part No.: 91196
Rapid clotting on body parts where tourniquets cannot be applied
Treats torso wounds

A hemostatic gauze pad is ideal to stop bleeding on areas of the body where you cannot apply a tourniquet like the torso, groin or neck. The pad is designed to promote rapid blood coagulation on an extreme bleeding wound. The pad is treated with an agent that absorbs the water in blood allowing the clotting factors to form a clot over the wound quickly.

Rapidly Distribute Critical Kits in Times of Crisis

This portable bag includes four specialized kits, designed to treat bleeding injuries on limbs, torso, and lungs in emergency situations when every second counts.



Contact your SSI Representative for all your First Aid Only Needs!

865-483-9332 or 800-229-7252 | customerservice@scisale.com



Respond to this email to unsubscribe.